


**PAJUNK®**

## *SonoSystem*

*The gold standard for  
outstanding echogenicity*



**Nerve blocks**


MADE IN GERMANY

## Cornerstone Reflectors

# Sono cannulas with maximum echogenicity

The visibility of cannula tips in ultrasound guided puncture is very important in preventing damage to neurons and blood vessels.<sup>1</sup> As even cannulas that are visible under ultrasound cannot always be identified at angles of 45° and above<sup>2</sup>, this property has become a key decision-making criterion in the selection of cannulas in practice.<sup>3</sup>


The patented Cornerstone Technology developed by PAJUNK® together with Dr. Chris Mitchell was designed specifically to solve this problem, and produces excellent visibility irrespective of the insertion angle.<sup>4</sup> Sono cannulas have a high degree of precision even at steep insertion angles. Both, the cannula shaft and -tip are very clearly visible.<sup>5</sup> In this way, Sono cannulas make an important contribution to the safety of the application.<sup>6</sup>



### Echogenic Cornerstone Geometry

The embossed structures in the Cornerstone Reflectors form three surfaces which meet each other at a 90° angle.


➔ This guarantees direct or indirect reflection of the ultrasound waves even at very steep insertion angles.<sup>2</sup>



### Sophisticated 360° arrangement

Both cannula segments are graduated all-around with evenly offset Cornerstone Reflectors. The number and layout of these reflectors is matched precisely to the relevant cannula diameter.

➔ Perfect cannula identification is guaranteed in every position.



### Visibility irrespective of the insertion angle

The Cornerstone Reflectors are designed that the ultrasound waves are very well reflected even with an insertion angle of 60° to 70°.<sup>2</sup>

➔ Ultrasound waves are reflected along a total length of 20 mm. Cannula shaft and cannula tip can be clearly identified.

1 Wiesmann et al., Compound imaging technology and echogenic needle ..., 2013; 38(5): 452–455




2 Uppal, Sondokoppam, Ganapathy, Effect of beam steering on ..., 2014; 61(10): 909–915

3 Sviggum, Ahn, Dilger, Smith, Needle echogenicity in sonographically ..., 2013; 32(1): 143–148

4 Hebard, Hocking, Echogenic technology can improve needle visibility ..., 2011; 36(2): 185–189

5 Edgcombe, Hocking, Sonographic identification of needle tip ..., 2010; 35(2): 207–211

6 Hocking, Mitchell, Optimizing the safety and practice ..., 2012, 604



- ➔ Echogenic, three-dimensional Cornerstone Geometry
- ➔ 360° configuration, arranged evenly around the cannula shaft
- ➔ Aligned to the cannula diameter
- ➔ Two 1 cm segments for positioning
- ➔ Cornerstone Reflectors are embossed as far as the tip of the cannula

= Optimum cannula visibility from shaft to tip, irrespective of the insertion angle



- ➔ Special configuration of the Cornerstone Reflectors for optimum 360° sonographic visibility
- ➔ Ultrasound waves are reflected along a length of 20 mm
- ➔ Clear identification of cannula shaft and tip
- ➔ Reflexion especially at steep insertion angles
- ➔ 360° length graduation for optimum positioning

## Dual guidance

# Precise nerve stimulation with NanoLine

The echogenic Cornerstone Cannula is used to identify the relevant nerve under ultrasound, and in a second step the "dual guidance" procedure can be performed to check the accuracy of catheter placement by means of nerve stimulation. To this end, the distance from the nerve is deduced from the lowest current strength required for stimulation. The NanoLine NanoLine thin layer technology developed by PAJUNK® offers considerable benefits with regard to accuracy of stimulation, because it allows the insulating layer to be reduced to a minimum without affecting the functionality. This extremely thin polymer layer, which is applied to every internal and external part of the device except the bare tip, allows highly accurate puncture and stimulation.

- ➔ Combination of ultrasound and stimulation technique
- ➔ Optimisation of puncture accuracy
- ➔ Better safety in use
- ➔ Nerve stimulator MultiStim ECO, designed specially for combination procedures
- ➔ Accurate stimulation and excellent gliding properties with NanoLine (only from PAJUNK®)



### Echogenic cannula tip

The facet grinding has two inclination angles, and is coated with NanoLine technology, apart from the tip, which is bare.

- ➔ Optimum conditions for outstanding cannula tip visibility.



### Precise stimulation

The very thin NanoLine coating guarantees a 100% insulation. The contact point at the cannula tip remains bare.

- ➔ Stimulation takes place only via the electroconductive puncture tip, generating a highly precise electrical field.



MultiStim ECO is a compact nerve stimulator developed by PAJUNK®, an easy-to-use device that meets the demands of combined procedures.



### Coated inner lumen

The thin-coating technology used, means that even the inner lumina of cannulas can be coated.

➔ This smooths out any unevenness and allows better flow of the anaesthetic.

### Reduced puncture force – increased glide properties

The very thin coating means that the exterior diameter is no different from when conventional coating techniques are used. It also produces an extreme surface smoothness.

➔ NanoLine cannulas glide easily through tissue and do not require great puncture force.

### The advantages of NanoLine:

- ➔ Layer thickness is reduced to a minimum
- ➔ There is no change to the external diameter of the cannula
- ➔ The same excellent insulation properties as with conventional procedures
- ➔ A smooth surface to reduce the puncture force
- ➔ Extremely accurate stimulation via the contact point at the tip of the cannula

## The kit variants

# Continuous peripheral nerve block anaesthesia under ultrasound

Together with Dr Meier, PAJUNK® is the first manufacturer to have developed a technique by which the catheter is introduced in sterile condition directly from a container through the cannula, and the anaesthetic is injected through the catheter. This patented technique has established itself successfully on the market and has aroused a great response and recognition in the professional world.

PAJUNK® differentiates between four SonoSystem kits depending on the catheter properties:

### SonoLong Echo for convincing sonographic visibility

Kit consisting of SonoLong Echo catheter + SonoLong NanoLine cannula  
optional with three different tip designs



Facet grinding



SPROTTE® SPECIAL tip




Tuohy tip




### SonoLong Curl Echo for precise catheter placement

Kit consisting of SonoLong Curl Echo catheter with coiled tip + SonoLong NanoLine cannula  
with Tuohy tip



Tuohy tip



You also have a choice between up to three different cannula types depending on the kit variant. While the facet grinding cannula exclusively places the catheter parallel to the nerve, the SPROTTE® SPECIAL cannula and the Tuohy cannula are suitable for those cases, where it is necessary to introduce the catheter at an angle to the nerve.

### SonoLong Sono for maximum kink resistance

Kit consisting of SonoLong Sono catheter with integrated stainless steel helical coil + SonoLong NanoLine cannula optional with three different tip designs



Facet grinding



SPROTTE® SPECIAL tip



Tuohy tip



### StimuLong Sono II for additional nerve stimulation

Kit consisting of StimuLong Sono catheter with stimulateable tip + SonoLong NanoLine cannula with facet grinding or Tuohy tip



Facet grinding



Tuohy tip



## SonoLong Echo

# Catheter placement directly through the cannula

The ultrasound procedure has radically changed processes in continuous regional anaesthesia. This has resulted in new requirements for the visibility of cannulas and catheters as well as their handling that our development department has been intensively working on for some time. Since PAJUNK® set a milestone in the industry in terms of the sonographic visibility of cannulas with Cornerstone Technology, we are pleased that with SonoLong Echo we succeeded in developing an ultrasound visible, MRT visible and radiopaque catheter.

### The kit consists of:

- ➔ SonoLong NanoLine cannula with Cornerstone Reflectors optionally available with SPROTTE® SPECIAL tip, Tuohy tip or facet grinding
- ➔ SonoLong Echo catheter with catheter container
- ➔ FixoLong system
- ➔ Colour coded adapter



### Ascending length indication

On the 50 cm long catheter an ascending length indication up to a length of 30 cm in intervals of 5 cm is provided. Its exact position can therefore be determined at any time.

### Simple handling

The SonoLong Echo kits are equipped with a catheter container. The catheter can therefore be introduced into the cannula in a sterile condition. Furthermore this prevents the common memory effect that arises when winding the catheter in the packaging.





The special catheter material is extremely visible under ultrasound monitoring and in MRT and also radiopaque. – The best requirements for clear identification in all three procedures.



### Steel stylet

The catheter of the SonoLong Echo kit has a steel stylet, which is locked in place in the introductory aid at its end, and is removed together with the container after the application of the catheter.

➔ This provides the catheter with outstanding stability.



### Catheter with central opening

The centre of the catheter is open.

➔ This permits the free flow of the anaesthetic – particularly in connection with the post-operative injection pump.

➔ *Under ultrasound monitoring and MRT clearly visible catheter material that is also radiopaque*

➔ *Is introduced into the cannula in a sterile condition*

➔ *Central opening enables a good and continuous flow of anaesthesia*


## SonoLong Curl Echo

# The “on-the-dot” anaesthesia, with a minimum of medication required

The provision of Cornerstone Reflectors eases precise positioning of the SonoLong NanoLine cannula in the direct proximity of the nerve under ultrasound monitoring. Since the catheter will follow the path with the least resistance, and this does not always coincide with the neural structures, it is also necessary to monitor the position of the catheter exactly. Against this background, the SonoLong Curl Echo kit was developed by PAJUNK® together with Dr. Cedric Luyet, which has a catheter that has been designed in a very special manner. The SonoLong Curl Echo catheter enables precise catheter positioning for a minimum of required medication. As soon as the catheter passes through the opening of the cannula, the soft tip of the Curl catheter tip will roll up and therefore access the point where the cannula tip is positioned. This permits an extremely precise anaesthesia, with a minimum of medication required. The SonoLong Curl Echo catheter is visible under ultrasound and is also radiopaque.

### The kit consists of:

- SonoLong NanoLine cannula with Cornerstone Reflectors and Tuohy tip
- SonoLong Curl Echo catheter in the catheter container
- FixoLong system
- Colour coded adapter





The radiopaque SonoLong Curl Echo catheter is extremely visible under ultrasound. It has a curled end, a closed tip and six lateral openings. It can therefore be positioned extremely precisely and also provides an even distribution of anaesthetic.

For position control after catheter placement the catheter has been equipped with a 10 mm marking. Therefore it can be read off, whether the catheter is still in the correct position

Marking for position al control

SonoLong Curl Echo catheter


Exact catheter positioning

The SonoLong Curl Echo catheter curls in the direct proximity of the nerve. This permits a very precise administration of local anaesthetic at the intended location.

Ø 10 mm

SonoLong Curl Echo catheter with SonoLong NanoLine cannula with Tuohy tip

- ➔ *Curled tip enables precise placement*
- ➔ *Six lateral eyes allow an even distribution of the anaesthesia*
- ➔ *Integrated stainless steel helical coil for maximum buckle protection and unhindered flow*
- ➔ *Cannula and catheter are extremely visible under ultrasound*



### Integrated stainless steel helical coil

The stainless steel helical coil confers kinking resistancy to the catheter.

➔ This enables the unhindered flow of the anaesthetic, also and particularly over a longer period of time.

## SonoLong Sono

# The kink resistant catheter with stainless steel helical coil

The SonoLong Sono kit is a joint development of PAJUNK® and Dr Meier. It differs from the SonoLong NanoLine kit only in catheter design. It is equipped with an integrated stainless steel helical coil, and is suitable for long-term treatment in pain therapy and plexus anaesthesia. The SonoLong Sono kit is, just as the SonoLong Echo kit, available with three different types of cannulas: with a SPROTTE® SPECIAL tip, cannula with facet grinding or with a Tuohy tip. All three cannula variants are equipped with the proven Cornerstone Reflectors for improved visibility under ultrasound monitoring. The catheter material itself is characterized by its excellent echogenicity.

### The kit consists of:


- ➔ SonoLong NanoLine cannula with Cornerstone Reflectors optionally available with SPROTTE® SPECIAL tip, Tuohy tip or cannula with facet grinding
- ➔ SonoLong Sono catheter in catheter container
- ➔ FixoLong system
- ➔ Colour coded adapter

#### Ascending length indication

On the 50 cm long catheter an ascending length indication up to a length of 30 cm in intervals of 5 cm is provided. Its exact position can therefore be determined at any time.

#### Simple handling


The SonoLong Sono kits are equipped with a catheter container. The catheter can therefore be introduced into the cannula in a sterile condition. Furthermore this prevents the common memory effect that arises when winding the catheter in the packaging.



#### Visibility and orientation

The SonoLong Sono catheter is manufactured from echogenic material and graduated along the first 30 cm.

- ➔ This simplifies positioning and identification under ultrasound.




Catheter interior

Its integrated stainless steel helical coil provides the SonoLong Sono catheter with the highest degree of mobility and simultaneous kinking resistancy – an important aspect for continuous applications.



SPROTTE® SPECIAL tip



### Unobstructed flow of the anaesthesia

The integrated stainless steel helical coil provides the catheter with kinking resistancy.

➔ The best requirements for the unobstructed flow of anaesthesia for continuous applications and the safe connection of the catheter to the injection pump.

### Soft, flexible catheter tip

The steel stylet ends 1 cm before reaching the catheter tip, which leaves the tip soft and flexible.

➔ This raises the flexibility during the insertion of the catheter and reduces the risk of unintentional injuries of vessels.

➔ *Maximum mobility and kinking resistancy thanks to the stainless steel helical coil*

➔ *Soft, flexible catheter tip prevents unintentional injuries*

➔ *Cannula and catheter are extremely visible under ultrasound monitoring*


## StimuLong Sono II

# The combination of ultrasound and nerve stimulation provides double safety

The StimuLong Sono II combines the advantages of a stimulateable catheter with those of echogenic Cornerstone Cannulas in one kit. The StimuLong Sono II catheter is characterised by its stimulateable tip. Stimulation is achieved here using an additional electric conductor that remains in the catheter over the entire application period and thus also enables subsequent position control. The relevant nerve is first identified by the echogenic Cornerstone Cannula. In a second stage, the accuracy of the catheter placement can then be checked using nerve stimulation. The distance to the nerve can be derived at the minimum current strength that is required for stimulation. As the inner lumen is completely NanoLine coated when compared with conventional cannulas, comprehensive insulation can be assumed. Stimulation takes place exactly when the catheter tip comes out of the cannula. A secondary position control of the catheter tip is also possible for intermittent post-operative pain therapy.

### The kit consists of:

- ➔ SonoLong NanoLine cannula with Cornerstone Reflectors optionally available with Tuohy tip or cannula with facet grinding
- ➔ StimuLong Sono catheter in catheter container
- ➔ FixoLong system
- ➔ Colour coded adapter




### Graduation

On the 50 cm long catheter an ascending longitudinal graduation up to a length of 30 cm in intervals of 5 cm is provided. Its exact position can therefore be determined at any time.

### Simple handling


The SonoLong Sono kits are equipped with a catheter container. The catheter can therefore be introduced into the cannula in a sterile condition. Furthermore this prevents the common memory effect that arises when winding the catheter in the packaging.



The atraumatically rounded tip of the StimuLong Sono catheter is gold-plated. As a result, the highest possible conductivity and an excellent stimulation capability are guaranteed.



Tuohy tip



### Integrated stainless steel helical coil

The stainless steel helical coil provides the catheter with kinking resistancy.

➔ This ensures the unhindered flow of the anaesthetic, also and particularly over a longer period of time.

### Stimulateable catheter

The catheter has a continuous electrical circuit.

➔ This establishes the electrical connection between the ClampingAdapter and stimulateable catheter tip.

- ➔ *Reliable flow of anaesthesia thanks to integrated stainless steel helical coil*
- ➔ *Stimulateable catheter ensures double safety during placement*
- ➔ *Excellent conductivity and stimulation capability due to gold-plated catheter tip*

## Innovative catheter fixation

# FixoLong and FixoCath – ensure freedom of movement

PAJUNK® developed two different solutions for catheter fixation on patients; FixoLong and FixoCath for catheter sizes 19 G and 20 G. Especially in continuous application, they prevent an accidental removal of the catheter if the patient moves or an obstruction of the anaesthetic flow by an unfavourable position.

### FixoLong

With FixoLong, the catheter and filter are fixed close to the catheter exit, enabling greater freedom of movement in all continuous applications.

Art. No. 001151-40



### FixoCath

FixoCath is simultaneously a transparent adhesive and fixation and is placed directly on the exit point. This also ensures significant freedom of movement for the patient.

Art. No. 001151-37Z



### Filter 0.2µm

The 0.2 µm bacterial filter prevents the passage of bacteria.


Art. No. 001151-37X



### ClampingAdapter

The PAJUNK® ClampingAdapter is a special anchoring device which prevents overstretching and ensures optimum anaesthetic flow.


Art. No. 001151-37V





# All information at a glance

## SonoLong Echo




- Kit consisting of:
- SonoLong NanoLine cannula with Luer Lock connector and electrical connecting cable
  - Adaptable injection tube
  - SonoLong Echo catheter 20 G x 50 cm with central opening and steel stylet
  - ClampingAdapter (yellow)
  - Filter 0.2 µm
  - FixoLong

Product	Size	SonoLong Echo catheter	Item No.	PU	
<b>SonoLong NanoLine cannula</b>					
	Facet grinding	19 G x 50mm	20 G x 50cm	531185-31A	10
		19 G x 75mm	20 G x 50cm	561185-31A	10
		19 G x 100mm	20 G x 50cm	521185-31A	10
		19 G x 150mm	20 G x 50cm	511185-31A	10
	SPROTTE® SPECIAL tip	19 G x 60mm	20 G x 50cm	531185-31B	10
		19 G x 120mm	20 G x 50cm	521185-31B	10
	Tuohy tip	18 G x 50mm	20 G x 50cm	531185-31C	10
		18 G x 75mm	20 G x 50cm	561185-31C	10
		18 G x 100mm	20 G x 50cm	521185-31C	10
		18 G x 150mm	20 G x 50cm	511185-31C	10


## SonoLong Curl Echo



- Kit consisting of:
- SonoLong NanoLine cannula with Luer Lock connector and electrical connecting cable
  - Adaptable injection tube
  - SonoLong Curl Echo catheter 20 G with integrated stainless steel helical coil, closed tip and 6 lateral openings
  - ClampingAdapter (yellow)
  - Filter 0.2 µm
  - FixoLong

Product	Size	SonoLong Curl Echo catheter	Item No.	PU	
<b>SonoLong NanoLine cannula</b>					
	Tuohy tip	18 G x 50mm	20 G x 50cm	531188-31C	10
			20 G x 90cm	531188-34C	10
		18 G x 100mm	20 G x 90cm	521188-34C	10

## SonoLong Sono





- Kit consisting of:
- SonoLong NanoLine cannula with Luer Lock connector and electrical connecting cable
  - Adaptable injection tube
  - SonoLong Sono catheter 20 G with integrated stainless steel helical coil, steel stylet and central opening
  - ClampingAdapter (yellow)
  - Filter 0.2 µm
  - FixoLong

Product	Size	SonoLong Sono catheter	Item No.	PU	
<b>SonoLong NanoLine cannula</b>					
	Facet grinding	19 G x 50mm	20 G x 50cm	531187-31A	10
		19 G x 75mm	20 G x 50cm	561187-31A	10
		19 G x 100mm	20 G x 50cm	521187-31A	10
		19 G x 150mm	20 G x 50cm	511187-31A	10
	SPROTTE® SPECIAL tip	19 G x 60mm	20 G x 50cm	531187-31B	10
		19 G x 120mm	20 G x 50cm	521187-31B	10
	Tuohy tip	18 G x 50mm	20 G x 50cm	531187-31C	10
		18 G x 75mm	20 G x 50cm	561187-31C	10
		18 G x 100mm	20 G x 50cm	521187-31C	10
		18 G x 150mm	20 G x 50cm	511187-31C	10

## StimuLong Sono II



- Kit consisting of:
- SonoLong NanoLine cannula with Luer Lock connector and electrical connecting cable
  - Adaptable injection tube
  - StimuLong Sono catheter 20 G x 50cm with central opening, electrically conductive stylet and integrated metal helical coil
  - StimuLong ClampingAdapter (yellow) with integrated stimulation connection
  - Connecting cable
  - Filter 0.2 µm
  - FixoLong

Product	Size	StimuLong catheter	Item No.	PU	
<b>SonoLong NanoLine cannula</b>					
	Facet grinding	19 G x 50mm	20 G x 50cm	531187-32A	10
		19 G x 100mm	20 G x 50cm	521187-32A	10
	Tuohy tip	18 G x 50mm	20 G x 50cm	531187-32C	10
		18 G x 100mm	20 G x 50cm	521187-32C	10




# SonoCover



SonoCover

with adhesive strips  
with two elastic bands  
with ultrasound gel  
with 3D chamber  
sterile

Product	Size	Item No.					
SonoCover	15 x 30 cm	021151-1530	•	•	•		•
SonoCover	15 x 60 cm	021151-1560	•	•	•		•
SonoCover	15 x 100 cm	021151-1501	•	•	•		•
SonoCover	15 x 30 cm	011151-1530	•	•		•	•
SonoCover	15 x 60 cm	011151-1560	•	•		•	•
SonoCover	15 x 100 cm	011151-1501	•	•		•	•
SonoCover	15 x 30 cm	031151-1530	•	•	•	•	•
SonoCover	15 x 60 cm	031151-1560	•	•	•	•	•
SonoCover	15 x 100 cm	031151-1501	•	•	•	•	•
SonoCover	15 x 30 cm	001151-1530	•	•			•
SonoCover	15 x 60 cm	001151-1560	•	•			•
SonoCover	15 x 100 cm	001151-1501	•	•			•



www.raguru.com

# Studies

- **Abbal B., Choquet O., Gourari A., Bouic N., Massone A., Biboulet P., Bringuier S., Capdevila X.** Enhanced visual acuity with echogenic needles in ultrasound-guided axillary brachial plexus block, *Minerva Anesthesiol.* 2015 April; 81(4): 369–378
- **Bischoff J. M., Koscielniak-Nielsen Z. J., Kehlet H., Werner M. U.** Ultrasound-guided ilioinguinal/iliohypogastric nerve blocks for persistent inguinal postherniorrhaphy pain: a randomized, double-blind, placebo-controlled, crossover trial, *Anesth. Analg.* 2012 Jun; 114(6): 1323–1329
- **Fuzier R., Casalprim J., Bataille B., Harper I., Magues J. P.** The echogenicity of nerve blockade needles, *Anesth.* 2015; 70: 462–466
- **Edgcombe H., Hocking G.** Sonographic identification of needle tip by specialists and novices: a blinded comparison of 5 regional block needles in fresh human cadavers, *Reg. Anesth. Pain Med.* 2010 March–April; 35(2): 207–211
- **Hebard S., Hocking G.** Echogenic technology can improve needle visibility during ultrasound-guided regional anesthesia, *Reg. Anesth. Pain Med.* 2011 March–April; 36(2): 185–189
- **Hebard S., Hocking G., Murray K.** Two-dimensional mapping to assess direction and magnitude of needle tip error in ultrasound-guided regional anaesthesia, *Anaesth. Intensive Care* 2011; 39(6): 1076–1081
- **Hocking G. Mitchell C.** Optimizing the safety and practice of ultrasound-guided regional anesthesia: the role of echogenic technology, *Curr. Opin. Anaesthesiol.* 2012 Oct; 25(5): 603–609
- **Lurf M., Leixnering M.** Ultraschallgezielte Anlage eines Schmerzkatheters am N. medianus am Unterarm. Schmerzfreie Physiotherapie nach Arthrolyse und Tenolyse, *Der Anaesthesist*, Volume 57 Number 7, 2008; 686–688 [Plexolong]
- **Luyet C., Meyer C., Herrmann G., Hatch G. M., Ross S., Eichenberger U.** Placement of coiled catheters into the paravertebral space, *Anaesth.* 2012; 67: 250–255
- **Luyet C., Seiler R., Herrmann G., Hatch G. M., Ross S., Eichenberger U.** Newly Designed, Self-Coiling Catheters for Regional Anesthesia – An Imaging Study, *Reg. Anesth. Pain. Med.*, Volume 36 Number 2, 2011; 171–176
- **Mahmoud K. M., Ammar A. S.** Ultrasound-guided continuous infraclavicular brachial plexus block using bupivacaine alone or combined with adenosine for pain control in upper limb surgery, *Saudi Journal of Anaesth.*, Volume 5 Number 2, 2011; 132–137 [Plexolong]
- **Morath U., Luyet C., Spadavecchia C., Stoffel M. H., Hatch G. M.** Ultrasound-guided retrobulbar nerve block in horses: a cadaveric study, *Vet. Anaesth. Analg.* 2013; 40(2): 205–211
- **Schafhalter-Zoppoth I., McCulloch C. E., Gray A. T.** Ultrasound Visibility of Needles Used for Regional Nerve Block: An In Vitro Study, *Reg. Anesth. Pain. Med.*, Volume 29 Number 5, 2004 Sept–Oct; 480–488
- **Schummer W., Sakka S. G., Hüttemann E., Reinhart K., Schummer C.** Ultraschall und Lagekontrolle bei der Anlage zentraler Venenkatheter, *Anaesthesist* 2009; 58: 677–685 DOI 10.1007/s00101-009-1569-1
- **Sviggum H. P., Ahn K., Dilger J. A., Smith H. M.** Needle echogenicity in sonographically guided regional anesthesia: blinded comparison of 4 enhanced needles and validation of visual criteria for evaluation, *J. Ultrasound Med.* 2013 Jan; 32(1): 143–148
- **Taboada M., Rodríguez J., Amor M., Sabaté S., Alvarez J., Cortés J., Atanassoff P. G.** Is Ultrasound Guidance Superior to Conventional Nerve Stimulation for Coracoid Infraclavicular Brachial Plexus Block?, *Reg. Anesth. Pain. Med.*, Volume 34 Number 4, 2009; 357–360 [short beveled needle]
- **Tsui B. C. H., Tsui J.** Reusable phantom with feedback signal for ultrasound needle tip control, *Reg. Anesth. Pain Med.* 2011; 36(6): 630–631
- **Uppal V., Sondekoppam R. V., Ganapathy S.** Effect of beam steering on the visibility of echogenic and non-echogenic needles: a laboratory study, *Can. J. Anesth.* 2014 Oct; 61(10): 909–915
- **Wiesmann T., Borntträger A., Zoremba M., Neff M., Wulf H., Steinfeldt T.** Compound imaging technology and echogenic needle design: effects on needle visibility and tissue imaging, *Reg. Anesth. Pain Med.* 2013 Sep–Oct; 38(5): 452–455